

TABLE OPTIONS CATALOGUE

November 2024





Welcome to the new Shor-Line Ltd Animal Housing Product Catalogue November 2024

For over 39 years, Shor-Line Ltd has maintained the quality of our products and customer service in its business philosophy.

As world leading brand suppliers to Veterinary Practitioners, we constantly refine and add to an extensive array of products to help you, the Veterinary professional in your vital role providing the best possible care for your patients.

So, please enjoy our new catalogue of products and we look forward to working with you.

Your sincerely,

The **SHOR-LINE** team.



SURGERY TABLES

CONTINUUM FLAT TOP SURGERY TABLE



Continuum Flat Top Hydraulic Base Surgery Table

- Maximum working height of table top 1041mm
- Minimum working height of table top 775mm
- Maximum load weight 113kg
- All stainless steel top with fully rounded corners
- Slight centre elevation with shallow trough around the table top edges
- Near vertical positioning
- Attachment rail on all sides of the table
- Four sliding quick release cam cleat tie down attachments
- Table base is designed so that the surgeon and other surgical equipment can come into close proximity to the table
- Four levelling screws to allow you to adjust the table base




Product Code	Description	W (mm)	H (mm) min	H (mm) max	L (mm)
903.4000.04	Continuum Flat Top Hydraulic Base Surgery Table	559	755	1041	1524

CONTINUUM V-TOP SURGERY TABLE



Continuum V-Top Hydraulic Base Surgery Table

- All Stainless Steel top
- Maximum Height 1041mm
- Minimum Height 775mm
- Maximum load weight 113kg
- A 76.2mm wide stainless steel trough collects fluid runoff into a drain outlet
- top panel opens 270° for easy access to the trough for removal of surgical waste
- Near vertical positioning
- Attachment rail on either side of the table
- Four sliding quick release cam cleat tie down attachments
- Table base is designed so that the surgeon and other surgical equipment can come into close proximity to the table
- Four levelling screws to allow you to adjust the table base




Product Code	Description	W (mm)	H (mm) min	H (mm) max	L (mm)
903.4100.01	Continuum V-Top Hydraulic Base Surgery Table	482.6	775	1041	1473

CLASSIC FLAT TOP SURGERY TABLES



Classic Flat Top Surgery Table

- All Stainless Steel
- Maximum Height 996.9 mm
- Minimum Height 730.2 mm
- Maximum load weight 136kg
- Fully rounded corners on table top
- Slight centre elevation with shallow trough around the table top edges
- 360° revolving capability and near vertical positioning
- Attachment rail on all sides of the table
- Four sliding quick release cam cleat tie down attachments
- Available with foot operated Hydraulic or Electric Base
- Table Base has the option of offset wheels for ease of movement and excellent positioning



Surgery Table is 558.8 mm wide x 1524mm long x between 730.2 & 996.9 mm high


Product Code	Description
903.4000.01	Flat-Top Surgery Table with Hydraulic Lift base
903.4040.00	Flat-Top Surgery Table with Electric Lift Base

CLASSIC V-TOP SURGERY TABLES



Classic V-Top Surgery Table

- All Stainless Steel
- Maximum Height 1041.4mm
- Minimum Height 774.7mm
- Maximum load weight 136kg
- Two V-Top panels adjust independently from either end of the table
- 360° revolving capability and near vertical positioning
- Attachment rail on either side of the table
- Two sliding quick release cam cleat tie down attachments
- 762mm stainless steel trough centred under the two V-Top panels
- Available with foot operated Hydraulic or Electric Base
- Table Base has the option of offset wheels for ease of movement and excellent positioning



Classic V-Top Surgery Table

Surgery Table is 482.6 mm wide x 1473.2 mm long x between 774.7 & 1041.4 mm high

Product Code	Description
903.4100.02	V -Top Surgery Table with Hydraulic Lift base
903.4140.00	V-Top Surgery Table with Electric Lift Base

SURGERY TABLE ACCESSORIES

Wheels for Surgery Table

Featuring two castors, the Wheels for Surgery Table can be easily fitted to your existing Shor-Line Surgery Table.

- Surgery Table can be moved by lifting the opposite side to the wheels and pushing the table top
- Stability is not an issue as the Surgery Table remains on its levelling legs until the table is lifted



Wheels for Surgery Table

Product Code	Description
803.0002.02	Wheels for Surgery Tables

Anti-Fatigue Mats

The Anti-Fatigue Mat is an ideal aid for minimising leg and back discomfort from standing during long procedures or surgeries.

- Features a dual-density closed cell foam
- Help to distribute weight without making you feel "unstable"
- The mat is 12mm thick with tapered and rolled edges which and can easily be cleaned with soap and water



Anti-Fatigue Mats

Product Code	Description
803.0100.05	Anti Fatigue Mat 48" x 24" (Pad dimensions 1219mm x 609mm)
803.0100.07	Anti Fatigue Mat 36" x 24" (Pad dimensions 914mm x 609mm)


SURGERY TABLE ACCESSORIES

Dental Extension

The Dental Extension allows standard Flat-Top and V-Top tables to become a wet table for procedures such as dentals and also provides extra space when you need a surgery table to be longer.



- All parts are heavy-gauge Stainless Steel



Height Extension

Product Code

Description

803.4001.08

Flat-Top Table Dental Extension

803.4100.08

V-Top Table Dental Extension

Collecting Pail

This large 11.8L capacity Collecting Pail features a free swinging handle for use with all surgery tables designed to hold collecting pails.

- Made from 100% seamless, non-corrosive stainless steel

Collecting Pail

Product Code

Description

190.9000.04

Collecting Pail




ELECTRIC LIFT TABLES

Cojack-Plus Multi-Purpose Table

The Cojack-Plus Multi-Purpose Table meets all your patient lifting, exam, transport and surgical procedure needs. The footswitch operates a dual rechargeable 24 volt actuator system that lifts up to 113kg.

The Cojack-Plus Multi Purpose Table features:

- Stainless steel top with raised edges to prevent fluid run-off
- Drain opening with hook for a collection pail
- Two tie down restraint cleats on each side for secure animal transport
- Adjustable IV holder
- Four non-marking bumpers
- Four nylon swivel and locking castors provide easy manoeuvrability



Cojack-Plus Multi Purpose Table

Product Code	Description	W (mm)	H (mm) min	H (mm) max	L (mm)
903.3130.01	Cojack-Plus Multi Purpose Surgery and Procedure Table	558.8	355.6	1016.0	1524.00

Cojack-Plus Multi Purpose Table Sizes

Mobile Animal Lift Table


The Shor-Line Mobile Animal Lift Table is ideal for lifting injured and post operative patients, saving time and effort.



- Maximum Load 90.7kg
- Stainless Steel Top with raised edges to prevent fluid run-off
- Bronze load bearings ensure long life and quiet operation
- Four fixed tie-downs for animal restraint
- Four non-marking bumpers
- Full charge in just 240 minutes
- Four nylon swivel and locking castors provide easy manoeuvrability
- Available with scale


Mobile Animal Lift Table Sizes

Product Code	Description	W (mm)	H (mm) min	H (mm) max	L (mm)
903.3120.00	Mobile Animal Lift Table	508 (Overall 787)	273	950	1130 (overall 1320)
903.3120.02	Mobile Animal Lift Table with Scale				



EXAMINATION TABLES

FOLD-UP WALL MOUNT



When folded against the wall, the Fold-up Wall-Mount Examination Table only extends 165 mm leaving your room free for other purposes.

- Uses torsion spring technology which allows you to easily and safely raise and lower the table
- Adjustable counterbalance spring for fine adjustment
- Adjustable stop allows for exact levelling of the work top
- Stainless Steel one piece die-formed top has raised edges to prevent fluid run off
- Fluid retaining shelf collects fluid run off if table is lifted before it is cleaned
- Rounded corners let you work comfortably from any angle

Fold-up Wall-Mount Examination Table

Product Code	Description	W (mm)	L (mm)
903.1130.04	Fold-up Wall-Mount Examination Table Peninsula	558.8	1117.6
903.1130.05	Fold-up Wall-Mount Examination Table Lateral	558.8	1117.6



Lateral

Peninsula



 **midmark® FOLDING EXAM TABLE**







The unique fold-down design of our wall-mounted folding exam table features a stationary work surface perfect for charts, instruments and client education materials. Plus, a special gas-assist lifting mechanism makes raising and lowering the table a breeze.

Midmark tables are designed to enhance your workflow and make you, your patients and clients as comfortable as possible.

Optional Traction Mat sold separately.

- 1067mm L (406mm folded down) x 610mm W
- Patient weight capacity; 57kg



TUB TABLES

KNEE SPACE TUB TABLE



White Laminate Knee Space Tub Table

Dirty procedures in practice require extensive clean-up. All of these procedures become easier, cleaner and more efficient when performed on our Multi-Purpose Tub Table.

- Rounded knee-space which eliminates sharp corners and edges and improves comfort during procedures
- Recessed 2 piece grill makes removal and clean-up quicker and easier
- 3 mm PVC wrapped on all doors and drawer edges to prevent chipping and for additional moisture protection

White Laminate Knee Space Tub Table Sizes

Product Code	Description	W (mm)	H (mm)
904.4000.10	Knee Space Tub Table with Door & 3 Drawers 152.4mm deep Right Hand Knee Space 1540 mm long	634	922 1540
904.4000.11	Knee Space Tub Table with Door & 3 Drawers 152.4mm deep Left Hand Knee Space 1540 mm long	634	922 1540

Economy Deck Mounted Mixing Faucet


Featuring a wide radius gooseneck nozzle with aerator and automatic diverter, water is directed over the centre of the Tub Table where you need it most. Large 101.6mm long faucet handles allow activation with wrist or forearm.

Also included is a pre-rinse sprayer with 1727.2mm flexible stainless steel hose. A 50.8mm spray nozzle head directs water flow where you need it.



Economy Deck Mounted Mixing Faucet

Product Code	Description
804.0006.83	Economy Deck Mounted Mixing Faucet




Lavatory Faucet with Hose and Spray

The Lavatory Faucet with Hose and Spray features a 203.2mm swing spout and 1219.2mm vinyl hose with automatic spray diverter. Water is automatically diverted to the spray when handle is used.

Lavatory Faucet with Hose and Spray

Product Code	Description
804.0006.03	Lavatory Faucet with Hose and Spray

PREP & UTILITY TABLES



Mobile Instrument Table

This convenient, versatile table is sized to fit easily over all standard small animal surgery tables. The large work surface provides ample room for your instruments and/or monitoring equipment.

- Handle adjusts the table top height from 1016mm to 1346.2mm which neatly slides under the table when not in use
- Table work surface measures 1016mm long x 584.2mm wide
- Four castors are equipped with locking brakes which keep the table in the desired position

Mobile Instrument Table

Product Code	Description	L1 (mm)	L2 (mm)	W (mm)	H minimum (mm)	H maximum (mm)
903.2003.00	Mobile Instrument Table	1016.0	838.2	584.2	1016.00	1524.0

Mayo Stand

This Mayo Stand is compact in size but large enough for most surgical procedures.

- Legs and base are constructed of chrome plated steel tubing
- Tray is constructed of stainless steel.
- Height is adjustable from 863.6 mm to 1346.2 mm with easy slide extension posts.
- Two grey hard rubber castors allow for easy movement.



Mayo Stand

Product Code	Description
803.0010.00	Mayo Stand

TERMS AND CONDITIONS

1. All prices are quoted for Shor-Line product ex-works Cowbridge, exclusive of VAT.
2. Payment of account is due (14) days from the date of invoice.
3. Late payment of account will attract a surcharge, equivalent to 2% per month from the date that the invoice is due, (e.g., 14 days from the date of invoice).
4. All products remain the property of Shor-Line Ltd until paid for in full.

PRODUCT WARRANTIES:

All Shor-Line stainless steel products are warranted against manufacturer defects in materials and/or workmanship, whilst they are the property of the original purchaser - assuming proper care and maintenance has been carried out. All mechanical parts and/or electrical devices carry a one year warranty, unless specifically extended or limited in the product description, product inserts or operation manuals.

Some United States warranties, as stated in the operation manuals do not apply outside the United States. These are superseded by European warranties - please check for details.

GUARANTEE OF SATISFACTION:

If a customer is not fully satisfied with the quality of workmanship of a Shor-Line product, Shor-Line must be notified in writing within (10) days of receipt of goods.

Notification should be mailed to: Customer Services Department

Shor-Line Limited

Unit 36

Vale Business Park

Llandow

Cowbridge

Vale of Glamorgan

CF71 7PF

United Kingdom

Telephone: +44 (0) 1446 77 20 41

Website: www.shor-line.co.uk

E-mail: info@shor-line.co.uk

Proof and date of purchase will be required for your request to be carried out quickly and efficiently.

CUSTOMERS' RECEIPT OF GOODS

Customers are advised to check that the number of packages corresponds with the quantity stated on the "Carrier's Note". Any discrepancy should be recorded on the note.

Any damages or shortages must be advised to "Shor-Line Limited" within (2) days, and confirmed in writing within (7) days. No other claim will be considered.

Please note that the claim procedures that are stated on the exterior of the packing material only apply to the United States market.

PRODUCT RETURN AUTHORISATION:

Customers who need to return products must first obtain authorisation from "Shor-Line Customer Services" on 01446 772041. Products may be returned for the following reasons:

1. *The product is received in a damaged or sub-standard condition.*
2. *The product requires "service" or "repair"*

To assist in the return authorisation please provide the following information:

- a. *The name of the product to be returned and the product code number.*
- b. *The product serial number*
- c. *The reason for return*
- d. *The date of purchase and the Shor-Line Invoice Number*

Shor-Line's Customer Service Department will allocate an Authorisation Number, which will ensure that the repair, replacement and/or credit will be handled quickly and efficiently. Returns will be subject to a 10% re-stock fee.

NOTES

Index

Page Number

Surgery Tables

Continuum Flat Top Surgery Table.....	5
Continuum V-Top Surgery Table.....	6
Classic Flat Top Surgery Table.....	7
Classic V-Top Surgery Table	8
Surgery Table Accessories.....	9-10
Electric Lift Tables	11

Examination Tables

Fold-up Wall Mount.....	13
Midmark Folding Exam Table	14

Tub Tables

Knee Space Tub Table and Taps	16
-------------------------------------	----

Prep & Utility Tables

Mobile Instrument Table.....	18
Mayo Stand	18

Terms and Conditions	19
----------------------------	----

