

QUINTIGARE™

Plus

Disinfectant Cleaner

Test results



MICRO-ORGANISM	DILUTION	CONTACT TIME (MINS)	TESTED BY	TEST PROTOCOL	RESULT
VIRUCIDAL ACTIVITY					
Canine parvovirus	1:50	30	BluTest	EN 14476	PASS
Vaccinia virus (VR-1549)*	1:50	5	BluTest	EN 14675	PASS
Avian Flu (Virus H5N1)	1:10	30	VLA**	DEFRA	PASS
FUNGICIDAL/YEASTICIDAL ACTIVITY					
Aspergillus brasiliensis	1:50	30	Abbott Analytical	EN 1657	PASS
Candida albicans	1:50	30	Abbott Analytical	EN 1657	PASS
Trichophyton mentagrophtes	1:50	5	Abbott Analytical	EN 1657	PASS
Microsporum canis	1:50	5	Abbott Analytical	EN 1657	PASS
Trichophyton rubrum	1:100	30	Abbott Analytical	EN 1657	PASS
MYCOBACTERICIDAL ACTIVITY					
Mycobacterium marinum	1:100	30	Abbott Analytical	EN 14204	PASS
Mycobacterium bovis	1:10	30	Abbott Analytical	EN 14204	PASS
Mycobacterium avium	1:100	30	Abbott Analytical	EN 14204	PASS
BACTERICIDAL ACTIVITY					
Pseudomonas aeruginosa	1:50	30	Abbott Analytical	EN1656	PASS
Proteus vulgaris	1:50	30	Abbott Analytical	EN1656	PASS
Staphylococcus aureus	1:50	30	Abbott Analytical	EN1656	PASS
Enterococcus hirae	1:50	30	Abbott Analytical	EN1656	PASS
Campylobacter jejuni	1:50	30	Abbott Analytical	EN1656	PASS
Listeria monocytogenes	1:50	30	Abbott Analytical	EN1656	PASS
Streptococcus equi subsp equi	1:100	5	Abbott Analytical	EN1656	PASS
Aeromonas hydrophilia	1:50	30	Abbott Analytical	EN1656	PASS
Salmonella typhimurium	1:50	30	Abbott Analytical	EN1656	PASS
MRSA	1:50	30	Abbott Analytical	EN1656	PASS
Bordetella bronchiseptica	1:50	10	Abbott Analytical	EN1656	PASS
Yersinia enterocolitica	1:100	5	Abbott Analytical	EN1656	PASS
Providencia alcalifaciens	1:50	30	Abbott Analytical	EN1656	PASS
Escherichia coli	1:100	5	Abbott Analytical	EN1276	PASS

*This is a surrogate to cover ENVELOPED VIRUSES. These include Poxviridae, Herpesviridae, Filoviridae (Ebola, Marburg), Flavivirus, Hepatitis C Virus, Hepatitis Delta Virus, Influenza, Paramyxoviridae, Rubella Virus, Measles Virus, Rabies Virus, Coronavirus (SARS, MERS), Human Immunodeficiency Virus, Human T Cell Leukemia Virus, Hepatitis B virus **Veterinary Laboratories Agency is an executive agency of the Department for Environment, Food & Rural Affairs.

IMPORTANT: This information was correct at the time of publication. As Safe Solutions (the manufacturer) have a policy to continuously test their products it is important that you contact us if you are unsure of any of the information contained within this document or if you require a more recent test result.

SHOR-LINE
SINCE 1927

Quality Equipment...

An Investment for Life!